



## ProBasics K1 Wheelchair

### STANDARD WHEELCHAIR

The **ProBasics K1 Manual Wheelchair** combines style and durability at an affordable price. Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings. Supporting patient weights up to 250 pounds, the ProBasics K1 Standard Wheelchair is a popular choice among providers. Available with swing-away footrests or elevating legrests.



Flip-back armrests



Easy-to-operate wheel locks



Front riggings with heel loop

### FEATURES

- Slide tube steel frame in black powder coat finish
- Single axle positioning
- Easy-to-clean, comfortable vinyl upholstery resists mildew and bacteria
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops or elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- HCPCS\* Code: K0001

K0195 OR E0990 (elevating legrests only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"
Seat Depth:	16"
Seat-to-Floor Height <sup>1</sup> :	19.5"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24 ½"   26 ½"   28 ½"
Overall Depth:	31 ¼" (without riggings)
Weight Capacity:	250 lb
Product Weight <sup>2</sup> :	41 – 43 lb
Carton Dimensions:	34" x 12" x 41"
Shipping Weight:	45 – 50 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC11616DS	16" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070558	each
WC11816DS	18" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070565	each
WC12016DS	20" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070572	each
WC11616DE	16" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070589	each
WC11816DE	18" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070596	each
WC12016DE	20" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070602	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



## ProBasics K2 Wheelchair

### STANDARD HEMI WHEELCHAIR



The **ProBasics K2 Standard Hemi Wheelchair** is designed to meet a variety of patient profiles. With a patient weight capacity of 300 pounds, a dual axle for hemi-height adjustment and multi-position caster forks, it is one of the more versatile wheelchairs on the market today. Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings. Available with swing-away footrests or elevating legrests.



Flip-back armrests



Easy-to-operate wheel locks



Front riggings with heel loop

### FEATURES

- Slide tube, powder-coated silver vein steel frame
- Dual axle provides easy transition of seat height to hemi-level
- Easy-to-clean, comfortable vinyl upholstery resists mildew and bacteria
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops or elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters
- HCPCS\* Code: K0002, K0001

K0195 OR E0990 (elevating legrests only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"
Seat Depth:	16"
Seat-to-Floor Height <sup>1</sup> :	Adult 19.5"   Hemi 17.5"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24 ½"   26 ½"   28 ½"
Overall Depth:	31 ¼"
Weight Capacity:	300 lb
Product Weight <sup>2</sup> :	41 – 43 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	45 – 50 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC21616DS	16" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070633	each
WC21816DS	18" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070640	each
WC22016DS	20" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070657	each
WC21616DE	16" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070664	each
WC21816DE	18" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070671	each
WC22016DE	20" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070688	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



## ProBasics K3 Wheelchair

### LIGHTWEIGHT WHEELCHAIR

The **ProBasics K3 Lightweight Wheelchair** offers all the benefits of a standard K3 wheelchair with the added feature of removable flip back height adjustable arms and 2-inch upholstery seat extension. Weighing less than 36 pounds, the ProBasics K3 also comes standard with tool-free adjustable front riggings and easy push-lock wheel locks.



Flip-back armrests are height adjustable



2" upholstery seat extension



Front riggings with heel loop

### FEATURES

- Cross-brace, powder-coated silver vein steel frame
- Removable flip back height adjustable arms allow for easy transfer
- 2-inch removable seat extension
- Dual axle provides easy transition of seat height to hemi-level
- Easy-to-clean, comfortable nylon upholstery resists mildew and bacteria
- Includes tool-free swing-away footrests with heel loops or elevating legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters
- HCPCS\* Code: K0003  
K0195 OR E0990 (elevating legrests only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"
Seat Depth:	16" - 18"
Seat-to-Floor Height <sup>1</sup> :	Adult 20"   Hemi 17"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24"   26"   28"
Overall Depth:	31" (without riggings)
Weight Capacity:	300 lb
Product Weight <sup>2</sup> :	35 - 38 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	40 - 46 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC31616DS	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070695	each
WC31816DS	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070701	each
WC32016DS	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070718	each
WC31616DE	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070725	each
WC31816DE	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070732	each
WC32016DE	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070749	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



# ProBasics K4 Wheelchair

## HIGH STRENGTH LIGHTWEIGHT WHEELCHAIR

The **ProBasics K4 Wheelchair** is a high-strength lightweight extremely versatile wheelchair. With a built-in seat and upholstery extension combined with all of its other unique features that come standard, the ProBasics K4 easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance rental product at a moderate price.



Flip-back armrests are height adjustable



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



Quick release rear axle

### FEATURES

- Quick-release axles allow for tool-free removal of rear wheels for fast hemi-height adjustment
- Built-in seat rail and upholstery extension easily adjusts seat depth from 16" to 18" without any additional upholstery
- Flip-back padded arm rests are height adjustable from 10" to 14" to accommodate a seat cushion
- Adjustable height back and adjustable angle back from 5° - 20°
- Comfortable nylon upholstery resists mildew and bacteria
- Chart pocket on backrest for additional convenience
- Includes tool-free adjustable front riggings with composite angle adjustable footplates for more relaxed foot positioning
- Elevating padded leg rests are equipped with a pneumatic piston for infinite elevating positions
- Adjustable anti-tippers with push-pin flip-up wheels
- Easy push-to-lock wheel locks
- HCPCS\* Code: K0004

K0195 or E0990 (elevating legrests only)

E0971 (anti-tippers)

K0040 (angle adjustment footplates)

E2201 (non-standard seat frame - 20" chair only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"
Seat Depth:	16" - 18"
Seat-to-Floor Height <sup>1</sup> :	Adult 19.5"   Hemi 17.5"
Back Height:	16" - 18"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5"   26.5"   28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb
Product Weight <sup>2</sup> :	32 - 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 - 39 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC41616DS	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070756	each
WC41816DS	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070763	each
WC42016DS	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070770	each
WC41616DE	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070787	each
WC41816DE	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070794	each
WC42016DE	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070800	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary 1/4-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



## ProBasics K4 Transformer

### TRANSFORMS FROM K4 TO A TRANSPORT CHAIR

The perfect upgrade from a standard K0004 wheelchair, the **ProBasics TRANSFORMER™** is a dual purpose, high strength lightweight wheelchair that quickly transforms from a self-propelled wheelchair to a compact transport chair with 8" wheels all around. Great as a transitional chair for accompanied excursions while retaining self-propel capabilities. Also allows the slimming down of the chair to navigate narrow doorways.



24-inch quick release rear wheels



Transforms into a Transport Chair



Flip back, desk length, padded arms

### FEATURES

- Lightweight aluminum frame
- 8-inch full swivel front casters mounted on aluminum forks
- 24-inch "quick release" rear wheels with composite handrims
- Bolted on 8-inch transport chair rear wheels
- Self-locking paddle release, flip-back padded desk length arms
- Lightweight black, woven nylon upholstery
- Fold down back with easy release paddles
- Complete with lightweight aluminum wheel locks for both 24-inch and 8-inch rear wheel modes
- Detachable flip-up rear anti-tippers
- Accessory pouch sewn into back
- Footrests include heel loops and composite footplates (optional elevating legrests sold separately - WCT4ELR)
- Folds for easy storage
- HCPCS\* Code: K0004

K0195 or E0990 (elevating legrests only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"
Seat Depth:	16"
Seat-to-Floor Height <sup>1</sup> :	Adult 20"   Hemi 17"
Back Height:	16"
Back Style:	Folds downs
Overall Width w/Arms:	18" - 29"
Overall Depth:	28" - 32" (without riggings)
Weight Capacity:	250 lb
Product Weight <sup>2</sup> :	32 - 34 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 - 37 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WCT41616DS	16" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070817	each
WCT41816DS	18" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070824	each
WCT42016DS	20" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070831	each
Item #	Description	UPC	UOM		
WCT4ELR	Optional Elevating Leg Rests	815067070848	pair		

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



## ProBasics Reclining Manual Wheelchair



The **ProBasics Reclining Manual Wheelchair** provides full back and headrest support as patient reclines from 90 to 160 degrees. The 13-inch headrest and padded, full-length arms are removable. Black, padded vinyl upholstery resists mildew and bacteria. The frame has a chip-resistant silver vein finish. Mag-style wheels are lightweight and set back on frame to prevent patient from tipping when in the reclined position.



Reclines from 90° – 160°



Caregiver activated hand brakes



13-inch headrest supports patient when reclined and is removable

### FEATURES

- Slide tube steel frame with an attractive silver vein finish that is chip resistant and easy to maintain
- Padded Desk-length, detachable armrests provide added patient comfort
- Easy-to-clean, comfortable vinyl upholstery resists mildew and bacteria
- Chart pocket on backrest provides additional convenience
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite positions and composite footplates
- Composite handrims with mag-style wheels are lightweight and maintenance free
- Wheels are set back on frame to prevent tipping
- Rear anti-tippers come standard
- Headrest extension comes standard with head immobilizer
- 8" x 1" solid rubber front casters
- Easy push-to-lock wheel locks
- HCPCS\* Code: K0001/K0006
  - E1226 (reclining back)
  - E0966 (headrest extension)
  - E0990 (elevating legrests only)

### SPECIFICATIONS

Seat Width:	16"   18"   20"   22"
Seat Depth:	17"
Seat-to-Floor Height <sup>1</sup> :	Adult 20"
Back Height:	19.5"   34"
Back Style:	Reclining
Overall Width w/Arms:	24½"   26½"   28½"
Overall Depth:	31¼" (without riggings)
Weight Capacity:	300 lb   400 lb (WCR2218E)
Product Weight <sup>2</sup> :	51 – 60 lb
Carton Dimensions:	40" x 12" x 36"
Shipping Weight:	71 – 75 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	Weight Limit	UPC	UOM
WCR1616E	16" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070893	Each
WCR1816E	18" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070909	Each
WCR2018E	20" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070916	Each
WCR2218E	22" x 17"	Removable, Desk-Length	Elevating Legrests	400 lb	815067070923	Each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



## ProBasics K7

### HEAVY DUTY WHEELCHAIR

The **ProBasics K7 Heavy Duty Wheelchair** features a robust design to accommodate patients with wider seat and higher weight capacity requirements. With seat widths from 24 to 28 inches, a weight capacity of 450 - 600 pounds and a heavy-duty slide tube frame, the ProBasics K7 offers a rugged yet attractive option for the bariatric patient population.



Removable, desk length, padded arms



Heavy duty vinyl upholstery



Aluminum footplate

#### FEATURES

- Heavy Duty, slide tube silver vein steel frame
- 8" x 2" semi-pneumatic casters are adjustable in 3 positions
- Dual axle provides easy transition of seat height to hemi-level
- 24" rear urethane wheels with durable hand rims
- Removable desk length arms with padded armrests
- Heavy duty vinyl padded upholstery with chart pocket
- Tool-free footrests include heel loops and Aluminum footplates
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite elevating positions and aluminum footplates
- Easy push-to-lock wheel locks
- HCPCS\* Code: K0007  
K00195 or E0990 (elevating legrests only)

#### SPECIFICATIONS

Seat Width:	22"   24"   26"   28"
Seat Depth:	18"   20"
Seat-to-Floor Height <sup>1</sup> :	Adult 21"   Hemi 19"
Back Height:	16"   18"
Back Style:	Fixed
Overall Width w/Arms:	30"   32"   34½"   37"
Overall Depth:	33"   34" (without riggings)
Weight Capacity:	450 - 600 lb
Product Weight <sup>2</sup> :	51 - 64 lb
Carton Dimensions:	35" x 13" x 38"
Shipping Weight:	58 - 70 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	Weight Limit	UPC	UOM
WC72218DS	22" x 18"	Removable, Desk-Length	Swing-Away Footrests	450 lb	815067070855	each
WC72418DS	24" x 18"	Removable, Desk-Length	Swing-Away Footrests	450 lb	815067070862	each
WC72620DS	26" x 20"	Removable, Desk-Length	Swing-Away Footrests	500 lb	815067072057	each
WC72820DS	28" x 20"	Removable, Desk-Length	Swing-Away Footrests	600 lb	815067072064	each
WC72218DE	22" x 18"	Removable, Desk-Length	Elevating Legrests	450 lb	815067070879	each
WC72418DE	24" x 18"	Removable, Desk-Length	Elevating Legrests	450 lb	815067070886	each
WC72620DE	26" x 20"	Removable, Desk-Length	Elevating Legrests	500 lb	815067072071	each
WC72820DE	28" x 20"	Removable, Desk-Length	Elevating Legrests	600 lb	815067072088	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

<sup>1</sup>NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

<sup>2</sup>NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.