

Steel Transport Chair

The **ProBasics Steel Transport Chair** is a lightweight steel-construction design with a chip-proof and easy-to-maintain finish. The flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport. Comes complete with seat belt, tool-free swing-away composite footrests and heel loops. Fits patient weights up to 300 pounds.











FEATURES

- Lightweight steel construction (only 25 lbs)
- Flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Stylish, silver vein carbon steel frame with permanent rear axel for strength and durability
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Nylon seat straps allow for ease of folding and lifting
- Push-to-lock rear wheel locks
- Composite 8"wheels with solid rubber tires are lightweight and maintenance free
- Tool-free, swing-away composite footrests with heel loops included
- HCPCS Code: E1038*

SPECIFICATIONS

51 2011 10/1110115			
Overall Width:	21.25" (open) 9" (closed)		
Width between Grips:	19"		
Overall Height:	37 "		
Overall Depth:	38.75" with riggings 26.75" without riggings		
Seat-to-Floor Height:	19"		
Back Height:	16"		
Seat Dimensions:	19" (W) x 16" (D)		
Patient Weight Capacity:	300 lb		
Product Weight:	25 lb		
Carton Dimensions:	23" x 11" x 31"		
Shipping Weight:	34 lb		
Warranty:	Limited lifetime on frame One year on all parts		

Item #	Description	UPC	Seat Width	Color	Wt. Cap.	UOM
TCS1916SV	ProBasics Steel Transport Chair	815067070480	19"	Silver Vein	300 lb	1/cs

^{*} Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.